

- Providing several visual target and visual function icons
- Small volume
- Easy to operate



## ACP-60

www.supore.com

### Liquid Crystal Eye Chart

- More clear images**
- Switching visual targets smoothly and quickly**
- Space-saving wall-hanging design**
- Thinner and more aesthetic design**

The input power supply parameters:

AC100-240V/50HZ/60HZ

Standard power consumption: About 45W

Power consumption: Less than 2W

Maximum screen resolution: 1440 × 900

Screen dot spacing: 0.284mm

Screen type: TFT

Screen ratio: 16:9 or 16:10 (slight difference is acceptable according to the actual batch)

Screen size: 21.5

Maximum brightness: 300cd/m

Maximum contrast: 2000:1

Signal response time: 5ms

Boundary dimension: 480 × 320 × 50mm

Use Distance: 1-7 m

Visual acuity record method:

The decimal English scores recording

Remote-controlled distance: No less than 8 m

Installation method: Bracket can be installed and hung on the wall

- Providing several visual target and visual function icons
- Small volume
- Easy to operate



## ACP-60S

www.supore.com

### Liquid Crystal Eye Chart

- More clear images**
- Switching visual targets smoothly and quickly**
- Space-saving wall-hanging design**
- Thinner and more aesthetic design**

The input power supply parameters:

AC100-240V/50HZ/60HZ

Standard power consumption: About 45W

Power consumption: Less than 2W

Maximum screen resolution: 1920\*1080

Screen dot spacing: 0.248mm

Screen type: TFT

Screen ratio: 16:9 or 16:10 (slight difference is acceptable according to the actual batch)

Screen size: 21.5

Maximum brightness: 300cd/m

Maximum contrast: 2000:1

Signal response time: 5ms

Boundary dimension: 515mm\*320mm\*50mm

Use Distance: 1-7 m

Visual acuity record method:

The decimal English scores recording

Remote-controlled distance: No less than 8 m

Installation method: Bracket can be installed and hung on the wall



# LX-88

## SLIT LAMP

www.supore.com

Warm LED light source, suitable for observing eye tissues and lesions  
"Frequency conversion" light control, light source color temperature is stable, always in the most suitable warm yellow light  
Humanized dimming mode, more convenient operation

Type : Slit Lamp  
Magnification : 6X 10X 16X 25X 40X  
Light Source : Warm LED

Slit Width : Continuous from 0 to 14mm  
Slit Length : Continuous from 1 to 14mm