

- Providing several visual target and visual function icons
- Small volume
- Easy to operate



ACP-60S

www.supore.com

Liquid Crystal Eye Chart

- More clear images**
- Switching visual targets smoothly and quickly**
- Space-saving wall-hanging design**
- Thinner and more aesthetic design**

The input power supply parameters:

AC100-240V/50HZ/60HZ

Standard power consumption: About 45W

Power consumption: Less than 2W

Maximum screen resolution: 1920*1080

Screen dot spacing: 0.248mm

Screen type: TFT

Screen ratio: 16:9 or 16:10 (slight difference is acceptable according to the actual batch)

Screen size: 21.5

Maximum brightness: 300cd/m

Maximum contrast: 2000:1

Signal response time: 5ms

Boundary dimension: 515mm*320mm*50mm

Use Distance: 1-7 m

Visual acuity record method:

The decimal English scores recording

Remote-controlled distance: No less than 8 m

Installation method: Bracket can be installed and hung on the wall