



- Adopting the latest national liquid crystal eye chart standard, containing 31 international visual targets including "E" letters, "C" letters and figures
- Stable operation
- Easy to operate



ACP-7

www.supore.com

Liquid Crystal Eye Chart

00000 0.5 ШЭШМ 0.5 КСАЕН 0.5 ✕✕✕✕ 0.6 КНШИБ 0.2 82635 0.5
 00000 0.6 ЕМЭШ 0.6 ZTOKV 0.6 ✕✕✕✕ 0.8 НМБШ 0.3 26389 0.6
 00000 0.8 МЕШЭ 0.8 CFLDN 0.7 ✕✕✕✕ 1.0 БНМК 0.4 89523 0.7

- Sophisticated appearance and delicate workmanship
- More accurate measurement according to the latest national liquid crystal eye chart standard
- 0.3 second response rate
- Various Kinds of Charts (A,B,C,D,F)

Icon:33selections

Speed of choosing icon:0.3Seconds
Display:1 exudes 3 horizontal lines
5 vertical lines 15 single letters

Filter color:Red/Green

Projector Distance:2.0m-7.0m(5m is Star

Image expansion of multiple:30×±5X (5m)

Projector size:252mm and 330(W)x225(H)mm Auto-standby function:3 minutes

Tilt:Ball joint

Light:6V 20W or (LED light), 2000 hours life

Power:AC220V 50Hz

Power:50VA

Size:Host:265(H)x230(W)x360(D)

Remote control:20(H)x61(W)x180(D)