

PL-3018II

Full Digital Portable Ultrasound Scanner




Features:

- 12.1 inch high resolution monitor
- More than 100 image permanent storage
- 8 segments digital TGC
- 2 probes sockets , support several optional probes
- Easy to operate and carry
- Measurement: distance, circumference, area, volume, heart rate, GA, FW, EDD.

Full Digital Portable Ultrasound Scanner

Create the health with science and technology



Specifications:

- Display depth (mm): 240 (max), 16 levels adjustable.
- Monitor size: 12 inch.
- Probe elements : 80.
- Display mode: B, B+B, B+M, M, 4B.
- Image gray scale: 256 level.
- Cine loop: 809 frame (max).
- Image storage: 32 frames.
- Scan angle: adjustable.
- Scan depth: 40mm-240mm.
- Acoustic power: 2 steps.
- Dynamic range: 100dB-130dB.
- Image flip: up/down, left/right, black/white.
- Focus position: adjustable.


- Focal space: 5 level.
- Measurement: distance, circumference, area, volume, heart rate, GA, FW, EDD.
- Notation: date, time, name, sex, age, doctor, hospital name. Full screen words edit.
- Output report: 2 type.
- Posture mark: ?40.
- Video output: PAL-D, VGA.
- Power Supply: AV220V±22V, 50MHz±1MHz.
- Machine dimension (mm): 320*165*290

Standard Configuration:

- Scanner Body: 1pc
- 3.5MHz Convex Probe: 1pc
- Charger with a Set of wires: 1set
- Bottle of USG gel: 1pc
- User Manual: 1pc

Optional:

- 5.0MHz Micro-convex Probe
- 6.5MHz Trans-vaginal Probe
- 7.5MHz Linear Probe
- Video Printer
- Trolley



Xuzhou Palmary Electronics Co., Ltd



Addr: No.4 Building, Zhongkai City Light, West Huaihai Road, Xuzhou City, Jiangsu Province, China
Tel: 86-516-83892808 , Fax: 86-516-83892898
E-mail: ultrasoundscanner@palmary.cn
Website: www.palmaryultrasound.com